

SBio HDL CHOLESTEROL Ppt. SET

(PEG Precipitation Method)

(For invitro diagnostic use only)

REF	90302010
Pack Size	10 ml



Store at 2-8°C	Manufacturer	This way up	L1 Precipitating Reagent	PEG Precipitation PEG Precipitation Method
Use by (Last day of stated month)	Consult Instructions for use	LOT Batch Number	IVD <i>In vitro</i> Diagnostic Medical Device	
Date of Manufacture	REF Catalogue Number	S HDL Cholesterol Standard (25 mg/dl)	EC REP Authorised Representative in the European Community	

INTENDED USE

HDL Cholesterol Ppt. SET is used for the determination of HDL Cholesterol in serum or plasma.

PRINCIPLE OF THE TEST

When the serum is reacted with the Polyethylene Glycol contained in the precipitating reagent, all the VLDL and LDL are precipitated. The HDL remains in the supernatant and is then assayed as a sample for cholesterol using the Cholesterol (CHOD/PAP) reagent.

CLINICAL SIGNIFICANCE

Lipoproteins are the proteins, which mainly transport fats in the blood stream. They can be grouped into chylomicrons, very low density lipoproteins (VLDL), low density lipoproteins (LDL) and high density lipoproteins (HDL). Chylomicrons and VLDL transport mainly triglycerides, though VLDLs also transport some amount of cholesterol. LDL carries cholesterol to the peripheral tissues where it can be deposited and increase the risk of arteriosclerotic heart and peripheral vascular disease. Hence high levels of LDL are atherogenic. HDL transports cholesterol from the peripheral tissues to the liver for excretion, hence HDL has a protective effect. The measurement of total and HDL cholesterol and triglycerides provide valuable information for the risk assessment of coronary heart diseases.

PRESENTATION

L1 : Precipitating Reagent 10 ml
S : HDL Cholesterol Standard (25 mg/dl) 5 ml

COMPOSITION

PEG 6000; Non Reactive Stabilizers, Detergents and Preservatives.

STORAGE / STABILITY

Contents are stable at 2-8°C till the expiry mentioned on the labels.

SAMPLE REQUIRED

Serum, EDTA plasma.

REAGENT PREPARATION

Reagents are ready to use.

SAMPLE WASTE AND DISPOSAL

Do not reuse the reagent containers, bottles, caps or plugs due to the risks of contamination and the potential to compromise reagent performance.

This product requires the handling of human specimen. It is recommended that all human sourced material are considered potentially hazardous and are handled in accordance with the OSHA

standard on blood borne pathogens.

Appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.

Handle specimen, solid and liquid waste and test components in accordance with local regulations and NCCLS guidelines M29, or other published biohazard safety guidelines.

PROCEDURE

Wavelength / filter : 505 nm (Hg 546 nm) / Green
Temperature : 37°C / R.T.
Light path : 1 cm

MATERIALS REQUIRED BUT NOT PROVIDED

General laboratory instrumentation like Spectrophotometer/Analyzer, Thermostatic Cuvette holder, Cuvettes, Micropipettes, Test tubes, Waterbath, Stopwatch/Timer.

Precipitation of VLDL & LDL:

Pipette into a clean dry test tube:

Precipitating Reagent (L1)	0.1 ml
Sample	0.1 ml

Mix well and incubate at R.T. for 5 min. Centrifuge at 2500-3000 rpm to obtain a clear supernatant.

Cholesterol Assay:

Pipette into clean dry test tubes labeled as Blank (B), Standard (S) and Test (T):

Addition Sequence	B (ml)	S (ml)	T (ml)
Working Reagent	1.0	1.0	1.0
Distilled Water	0.05	--	--
HDL Cholesterol Standard (S)	--	0.05	--
Sample	--	--	0.05

Mix well and incubate at 37°C for 5 min. or at R.T. (25°C) for 15 min. Measure the absorbance of the Standard (Abs. S), and Test Sample (Abs. T) against Blank, within 60 minutes.

CALCULATIONS

$$\text{HDL Cholesterol in mg/dl} = \frac{\text{Abs. T}}{\text{Abs. S}} \times 25 \times 2$$

(Where 2 is the dilution factor due to the deproteinization step)

Calculation of LDL Cholesterol (mg/dl): (Freidewald's Formula)

$$= (\text{Total Cholesterol}) - \left(\frac{\text{Triglycerides}}{5} \right) - (\text{HDL Cholesterol})$$

Freidewald's Formula is reliable provided that:

- No chylomicrons are present i.e., it is a fasting sample.
- Triglyceride values are below 400 mg/dl.
- Type III hyperlipoproteinemia is absent.

QUALITY CONTROL

The following process is recommended for QC during the assay of HDL Cholesterol. *Define and establish acceptable range for your laboratory.

- Two levels of control (Normal and Abnormal) are to be run on a daily basis.
- If QC results fall outside acceptance criteria, recalibration may be necessary.
- Review QC results and run acceptance criteria following a change of reagent lot.

SPECIFIC PERFORMANCE CHARACTERISTICS

Linearity:

This procedure is linear upto 150 mg/dl of HDL Cholesterol. If values exceed this limit, dilute the serum with normal saline (NaCl 0.9%) and repeat the assay. Calculate the value using the proper dilution factor.

Limit of detection:

The limit of detection for HDL Cholesterol is 0.5 mg/dl.

Interferences:

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Precision:

Precision studies were performed with two controls using NCCLS protocol EP5-A. The results of the precision studies are shown below:

Sample	Within-run		Between-run		Total	
	Mean	CV%	Mean	CV%	Mean	CV%
Control 1	39.3	6.63	39.5	4.58	78.8	11.21
Control 2	25.92	7.58	26.28	4.86	52.2	12.44

Method comparison:

Comparative studies were done to compare our reagent with another commercial HDL Cholesterol Assay. No significant differences were observed. Details of the comparative studies are available on request.

REFERENCE RANGE

		Prognostically Favourable	Standard Risk Level	Risk Indicator
HDL Chol. (mg/dl)	Males	>55	35-55	<35
	Females	>65	45-65	<45
LDL Chol. (mg/dl)	Males	<150	150-190	>190
	Females			
Total Chol. (mg/dl)	Males	>3.8	3.8-5.9	<5.9
HDL Chol. (mg/dl)	Females	>3.1	3.1-4.6	<4.6

It is recommended that each laboratory establish its own normal range representing its patient population*.

NOTE

HDL Cholesterol are reported to be stable in serum for 7 days when stored at 2-8°C. The sample should preferably be of 12 to 14 hours fasting. The supernatant should be clear. If it is hazy or cloudy, the sample should be diluted 1+1 with normal saline (NaCl 0.9%) and the precipitation step should be repeated. (Results x 2). Anticoagulants such as fluoride, oxalates and hemolysed serums should not be used. The reagent may be used in several automated analyzers. Instructions are available on request. Standard is traceable to standard reference material (SRM) 909b. Do not use turbid, deteriorated or leaking reagents.

REFERENCES

- Trinder P., (1969) Ann. Clin. Biochem. 6:24.
- Allain C.C., et. Al., (1974) Clin. Chem. 20:470.
- Flegg H.M., (1972) Ann. Clin. Biochem. 10:79.
- Grillo F., et. al. (1981) Clin. Chem. 27 : 375.
- Demacker P.N.M. et. al. (1980) Clin. Chem. 26 : 1775.
- Assmann G. (1982) Lipid Metabolism and Atherosclerosis, Schattauer Verlag, Stuttgart.



Mfd. for:

Singapore **SB**
Biosciences PTE Ltd.

11 Yishun Street 51, # 04-23, The Criterion, Singapore 767971.



CMC Medical Devices & Drugs S.L.,
C/ Horacio Lengo No. 18, CP 29006, Malaga, Spain.